



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Protective Guard

**Product Code :**  
LBNY-0005-1820005



**Description :**

To provide maximum protection to the face from splinters, etc., Protective Guard is supplied with a 450 mm square acrylic shield. Protective Guard is 2 metres high by 1.2 metres wide, and is fitted with wire mesh mounted on a stable frame. This shield is 5 mm thick and hooks onto the wire mesh in any position. Protective Guard offers a high degree of protection to the operators from flying debris, etc when used along side the MAGNUS Test frame.

**Technical Specification :**

A 450 mm square transparent plastic shield to hook onto the wire mesh guard  
A free standing wire mesh guard 2 m high and 1.2 m wide.



**Laboratory Instrument India**