



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Twin Hydraulic Ram System

**Product Code :**  
LBNY-0005-1820007



**Description :**

Both Twin Hydraulic Ram System are operated using a hand pump, which is mounted on a mobile table. The digital display is calibrated in Psi. The digital indicator can be moved around the hydraulic table and tilted to give a better viewing angle. A lever on the end of the hand pump allows direction control of the ram between raising and lowering. A digital pressure indicator shows the ram pressure via a transducer. Twin Hydraulic Ram System Two 20 tonne dual direction rams are each fitted in wheeled carriages that roll along the flanges of the top members of the MAGNUS frame. There are quick acting clamps that lock the carriage in position.

**Technical Specification :**

Rams operate from a hand pump mounted on a mobile table, with digital pressure indicator and integral shelf  
Two 20 tonne double acting rams on moveable carriages each to fit the 300 kN Universal Testing Frame  
To be used in conjunction with the 300 kN Universal Testing Frame and its standard set of accessories.  
Two hoses per ram with self sealing couplings to be included.



**Laboratory Instrument India**

