



Email : sales@laboratoryinstrumentindia.com

Product Name :
LLIN, 180x160x150cm LxWxH Polyester

Product Code :
LBNY-0006-2050001



Description :

Rectangular mosquito net, polyester 100 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, dimensions 180 x 160 x 150 cm (length x width x height) +/-5%.

Technical Specification :

Material: (composition ISO1833) 100% Polyester
Fabric: (construction ISO 8388) Warp knitted
Denier: (ISO 2060, DUPRO)100 denier +/-5%
Without door, with finished edge
Weight: (ISO 3801) 30-40 g/m² +/-10%
Netting Mesh size: Min. 24 holes/cm²
Dimensional stability to washing: (ISO5077/ISO6330) 8A at 30°C +/-5%
Hanging: 6 loops at the upper corners and in the middle
Insecticide: WHO PQ qualified
Concentration of insecticide: WHO PQ qualified
Target level of concentration: WHO PQ qualified.
Bursting strength: (ISO13938-1 or ISO 13938-2): Min. 350-405 kPa for net material and seams (ISO 13938-1 or ISO 13938-2)
(measuring area: 7.3 cm²)
Fire safety: (16 - CFR part 1610-CS191-53) Class 1.



Laboratory Instrument India