



Email : sales@laboratoryinstrumentindia.com

Product Name :
LLIN, 190x180x150cm LxWxH Polyester

Product Code :
LBNY-0006-2050002



Description :

Rectangular mosquito net, polyester 100 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, dimensions 190 x 180 x 150 cm (length x width x height) +/-5%.

Technical Specification :

Material: 100% Polyester
Fabric: Warp knitted
Denier: 100 denier +/-5%
Without door, with finished edge
Bursting strength: Min. 350-405 kPa for net material and seams (measuring area: 7.3 cm²)
Fire safety: Class 1
Target level of concentration: WHO PQ qualified
Concentration of insecticide: WHO PQ qualified
Weight: 30-40 g/m² +/-10%
Netting Mesh size: Min. 24 holes/cm²
Dimensional stability to washing: 8A at 30°C +/-5%
Hanging: 6 loops at the upper corners and in the middle
Insecticide: WHO PQ qualified.

Laboratory Instrument India

