



Email : sales@laboratoryinstrumentindia.com

Product Name :

LLIN, 190x180x180cm LxWxH Polyester

Product Code :

LBNY-0006-2050003



Description :

Rectangular mosquito net, polyester 100 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, dimensions 190 x 180 x 180 cm (length x width x height) +/-5%.

Technical Specification :

Material: 100% Polyester
Fabric: Warp knitted, monofilament
Denier: 100 denier, +/-5%
Without door, with finished edge.
Bursting strength: Min. 350 kPa for net material and seams
Fire safety: Class 1
Insecticide: WHO PQ qualified
Hanging: 6 loops at the upper corners and in the middle
Weight: Min. 40 g/m² +/- 10%
Netting Mesh size: Min. 24 holes /cm²
Dimensional stability to washing: 8A at 30°C +/-5%
Concentration of insecticide: WHO PQ qualified
Target level of concentration: WHO PQ qualified.

Laboratory Instrument India

