



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

LLIN, alt.dimensions LxWxH Polyester

**Product Code :**

LBNY-0006-2050004



**Description :**

Rectangular mosquito net, polyester 100 DENIER, treated with WHO PQ qualified insecticide, color white/blue/green, alternative dimensions.

**Technical Specification :**

Material: 100% Polyester

Fabric: Warp knitted

Denier: 100 denier +/-5%

Without door, with finished edge

Bursting strength: Min. 250,350-405 kPa for net material and seams (measuring area: 7.3 cm<sup>2</sup>)

Fire safety: Class 1

Hanging: 6 loops at the upper corners and in the middle

Weight: 30-40 g/m<sup>2</sup> +/-10%

Netting Mesh size: Min. 24 holes/cm<sup>2</sup>

Dimensional stability to washing: 8A at 30°C +/-5%



**Laboratory Instrument India**

