



Email : sales@laboratoryinstrumentindia.com

Product Name :
LLIN, 180x160x150cm LxWxH Polyethylene

Product Code :
LBNY-0006-2050005



Description :

Rectangular mosquito net, polyethylene 120-150 DENIER, treated with WHO PQ qualified insecticide, color white/blue/green, dimensions 180 x 160 x 150 cm (length x width x height) +/-5%.

Technical Specification :

Material: 100% Polyethylene
Fabric: Warp knitted, monofilament
Without door, with finished edge
Denier: 120-150 denier +/-5%
Fire safety: Class 1
Hanging: Min. 6 loops at the upper corner and in the middle.
Netting Mesh size: Min. 56 holes/inch²
Dimensional stability to washing: 8A at 30°C +/-5%
Bursting strength: Min. 350 kPa for net material and seams (measuring area: 7.3 cm²).



Laboratory Instrument India

