



Email : sales@laboratoryinstrumentindia.com

Product Name :

LLIN, 180x160x150cm LxWxH Polyethylene

Product Code :

LBNY-0006-2050005



Description :

Rectangular mosquito net, polyethylene 120-150 DENIER, treated with WHO PQ qualified insecticide, color white/blue/green, dimensions 180 x 160 x 150 cm (length x width x height) +/-5%.

Technical Specification :

Material: 100% Polyethylene

Fabric: Warp knitted, monofilament

Without door, with finished edge

Denier: 120-150 denier +/-5%

Fire safety: Class 1

Hanging: Min. 6 loops at the upper corner and in the middle.

Netting Mesh size: Min. 56 holes/inch²

Dimensional stability to washing: 8A at 30°C +/-5%

Bursting strength: Min. 350 kPa for net material and seams (measuring area: 7.3 cm²).



Laboratory Instrument India

