



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Welding Table

**Product Code :**  
LBNY-0004-820006



**Description :**

Welding Table

**Technical Specification :**

Dimensions(mm) 1000\*900\*820  
Frame made from steel tubes 50\*50\*1.7  
Top of table from steel mesh sheets 6mm thickness  
Table parts is very well assembled and supported  
Each table is inside cabin  
All cabins is inside cabin connected by ducts to exit with im Peller.



**Laboratory Instrument India**