



Email : sales@laboratoryinstrumentindia.com

Product Name :

LLIN, 190x180x150cm LxWxH Polyethylene

Product Code :

LBNY-0006-2050006



Description :

Rectangular mosquito net, polyethylene 120-150 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, dimensions 190 x 180 x 150 cm (length x width x height) +/-5%.

Technical Specification :

Without door, with finished edge

Material: 100% Polyethylene

Denier: 120-150 denier +/-5%

Fabric: Warp knitted, monofilament

Fire safety: Class 1

Hanging: Min. 6 loops at the upper corner and in the middle

Netting Mesh size: Min. 56 holes/inch²

Dimensional stability to washing: 8A at 30°C +/-5%

Bursting strength: Min. 350 kPa for net material and seams (measuring area: 7.3 cm²)



Laboratory Instrument India

