



Email : sales@laboratoryinstrumentindia.com

Product Name :
LLIN, 190x180x180cm LxWxH Polyethylene

Product Code :
LBNY-0006-2050007



Description :

Rectangular mosquito net, polyethylene 120-150 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, dimensions 190 x 180 x 180 cm (length x width x height) +/-5%.

Technical Specification :

Material: 100% High Density Polyethylene
Denier: 120-150 denier, +/-5%
Fabric: Warp knitted, monofilament
Without door, with finished edge.
Bursting strength: Min. 350 kPa for net material and seams
Fire safety: Class 1
Hanging: 6 loops at the upper corners and in the middle.
Weight: Min. 33 g/m² +/- 10%
Netting Mesh size: Min. 17 holes /cm²
Dimensional stability to washing: 8A at 30°C +/-5%



Laboratory Instrument India

