



Email : sales@laboratoryinstrumentindia.com

Product Name :

LLIN, alt.dimensions LxWxH Polyethylene

Product Code :

LBNY-0006-2050008



Description :

Rectangular mosquito net, polyethylene 120-150 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, alternative dimensions.

Technical Specification :

Material: 100% Polyethylene
Fabric: Warp knitted, monofilament
Without door, with finished edge
Denier: 120-150 denier +/-5%
Bursting strength: Min. 350 kPa for net material and seams (measuring area: 7.3 cm²)
Fire safety: Class 1
Hanging: Min. 6 loops at the upper corner and in the middle.
Weight: Min. 38 g/m² +/-10%
Netting Mesh size: Min. 56 holes/inch²
Dimensional stability to washing: 8A at 30°C +/-5%



Laboratory Instrument India

