



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

LLIN, alt.dimensions LxWxH Polyethylene

**Product Code :**

LBNY-0006-2050008



**Description :**

Rectangular mosquito net, polyethylene 120-150 DENIER, treated with WHO PQ qualified insecticide, colour white/blue/green, alternative dimensions.

**Technical Specification :**

Material: 100% Polyethylene  
Fabric: Warp knitted, monofilament  
Without door, with finished edge  
Denier: 120-150 denier +/-5%  
Bursting strength: Min. 350 kPa for net material and seams (measuring area: 7.3 cm<sup>2</sup>)  
Fire safety: Class 1  
Hanging: Min. 6 loops at the upper corner and in the middle.  
Weight: Min. 38 g/m<sup>2</sup> +/-10%  
Netting Mesh size: Min. 56 holes/inch<sup>2</sup>  
Dimensional stability to washing: 8A at 30°C +/-5%



**Laboratory Instrument India**

