



Email : sales@laboratoryinstrumentindia.com

Product Name :
ELISA, reader, 8 channel

Product Code :
LBNY-0007-2150001



Description :

Used in reference and referral testing and high throughput screening. Used for many types of ELISA. For clinical laboratory applications.

Technical Specification :

Bench-top, open system, automatic
Spectral range 400-900nm
Absorbance range 0.001-4.000 OD
Panel assay for up to 8-panel assays in a single run
Single, dual & multi-wavelength endpoint measuring modes
User programmable measurement criteria for validation and interpretation of results (cut-off, transformation, & validation)
Up to 8 standards per tests; NC, PC, and QC.
Up to 200 user-definable assay protocols.
Up to 1000 patients, 10000 sample records.
Endpoint and kinetic assay modes
Accommodates all common 96-well ELISA plates
Automatic and manual plate mapping of assay
Power requirements 220V/50Hz, with voltage surge protection.
Eight measuring plus one reference photodetector
Built-in shaking (low, medium, high, and variable speeds)
On-board qualitative & quantitative data reduction software in English, French, and Spanish.



Laboratory Instrument India