



Email : sales@laboratoryinstrumentindia.com

Product Name :
ELISA, washer, 8 channel

Product Code :
LBNY-0007-2150002



Description :

A washer for ELISA microplates designed to ensure thorough washing of reagents between ELISA steps. This equipment is used for many types of ELISA for clinical laboratory applications.

Technical Specification :

Bench-top, open and automatic system
8-channel strip manifold
Wash parameters include: 16-character assay name, number of cycles, wash volume, flow rate and variable soak times, aspiration and dispensing
Crosswise aspiration/double aspiration of flat bottom micro-plates for reduced residual liquid
Built-in vacuum & pressure pump assembly
Bottles for waste, rinse and wash solution
Wash cycles: Between 1 to 10
Dispensing volumes from 25 to 3000ul
Accommodates flat, round, U or V-shaped bottom plates
75 user-definable protocols
Soak time 1-600 seconds
Power requirements 110-220V / 50Hz-60Hz, with voltage surge protection
Spill-over protection & electronics isolated from fluidics
Flip-out aerosol cover or similar.



Laboratory Instrument India