



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
ELISA, washer, 8 channel

**Product Code :**  
LBNY-0007-2150002



**Description :**

A washer for ELISA microplates designed to ensure thorough washing of reagents between ELISA steps. This equipment is used for many types of ELISA for clinical laboratory applications.

**Technical Specification :**

Bench-top, open and automatic system  
8-channel strip manifold  
Wash parameters include: 16-character assay name, number of cycles, wash volume, flow rate and variable soak times, aspiration and dispensing  
Crosswise aspiration/double aspiration of flat bottom micro-plates for reduced residual liquid  
Built-in vacuum & pressure pump assembly  
Bottles for waste, rinse and wash solution  
Wash cycles: Between 1 to 10  
Dispensing volumes from 25 to 3000ul  
Accommodates flat, round, U or V-shaped bottom plates  
75 user-definable protocols  
Soak time 1-600 seconds  
Power requirements 110-220V / 50Hz-60Hz, with voltage surge protection  
Spill-over protection & electronics isolated from fluidics  
Flip-out aerosol cover or similar.



**Laboratory Instrument India**