




Email : sales@laboratoryinstrumentindia.com

Product Name :
Cleaner, for HemoCue photometers/PAC-5

Product Code :
LBNY-0007-2250001



Description :

Used to clean the optical lens of a HemoCue photometer series by absorbing dirty, dust or blood stains.
Material: Sponge made of polyvinyl acetate, pre-moistened with a cleaning solution.
Cleaner made of a lint and fiber-free material, highly absorbent.
Moistened with a cleaning solution; a soaked sponge with specific shape, mounted on a plastic holder.
Pack of 5 cleaners, each cleaner individually packed in a welded aluminum foil.



Laboratory Instrument India