




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cleaner, for HemoCue photometers/PAC-5

**Product Code :**  
LBNY-0007-2250001



**Description :**

Used to clean the optical lens of a HemoCue photometer series by absorbing dirty, dust or blood stains.  
Material: Sponge made of polyvinyl acetate, pre-moistened with a cleaning solution.  
Cleaner made of a lint and fiber-free material, highly absorbent.  
Moistened with a cleaning solution; a soaked sponge with specific shape, mounted on a plastic holder.  
Pack of 5 cleaners, each cleaner individually packed in a welded aluminum foil.



**Laboratory Instrument India**