



Email : sales@laboratoryinstrumentindia.com

Product Name :

Microcuvette, for Glucose 201+/BOX-100

Product Code :

LBNY-0007-2250002



Description :

Single-use, sample collection and measuring microcuvette.
Capillary, venous and arterial whole blood can be used.
Collects 5 to 10uL capillary whole blood directly from the skin surface by capillarity.
Box of 100 microcuvettes , in plastic vial, with screw cap.
Inserted for direct reading into S0000590, Photometer,HemoCue glucose 201+/set.
Supplied with instruction manual in English, French, Spanish.

Technical Specification :

Microcuvette, for Glucose 201+/BOX-100



Laboratory Instrument India