



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Microcuvette, for Glucose 201+/BOX-100

**Product Code :**

LBNY-0007-2250002



**Description :**

Single-use, sample collection and measuring microcuvette.  
Capillary, venous and arterial whole blood can be used.  
Collects 5 to 10uL capillary whole blood directly from the skin surface by capillarity.  
Box of 100 microcuvettes , in plastic vial, with screw cap.  
Inserted for direct reading into S0000590, Photometer,HemoCue glucose 201+/set.  
Supplied with instruction manual in English, French, Spanish.

**Technical Specification :**

Microcuvette, for Glucose 201+/BOX-100



**Laboratory Instrument India**