



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Microcuvette, for Hb 301/BOX-200

**Product Code :**  
LBNY-0007-2250005



**Description :**

Capillary, venous and arterial whole blood can be used.  
Inserted for direct reading into S0000591, Photometer, HemoCue Hb 301/set.  
Microcuvette is made of polystyrene plastic, contains pre-applied, dried reagents.  
Single-use, sample collection and measuring microcuvette.  
Collects 10 to 12uL capillary whole blood directly from the skin surface by capillarity  
Box of 200 microcuvettes, 4x50, in a plastic vial, with screw cap.  
Supplied with instruction manual in English, French, Spanish.  
Supplier's product reference: HemoCue Hb 301 Microcuvettes.

**Technical Specification :**

Microcuvette, for Hb 301, box/200



**Laboratory Instrument India**