



Email : sales@laboratoryinstrumentindia.com

Product Name :
Microcuvette, for Hb 301/BOX-200

Product Code :
LBNY-0007-2250005



Description :

Capillary, venous and arterial whole blood can be used.
Inserted for direct reading into S0000591, Photometer, HemoCue Hb 301/set.
Microcuvette is made of polystyrene plastic, contains pre-applied, dried reagents.
Single-use, sample collection and measuring microcuvette.
Collects 10 to 12uL capillary whole blood directly from the skin surface by capillarity
Box of 200 microcuvettes, 4x50, in a plastic vial, with screw cap.
Supplied with instruction manual in English, French, Spanish.
Supplier's product reference: HemoCue Hb 301 Microcuvettes.

Technical Specification :

Microcuvette, for Hb 301, box/200



Laboratory Instrument India