



Email : sales@laboratoryinstrumentindia.com

Product Name :
Marker, diamond, set/3

Product Code :
LBNY-0007-2260001



Description :

Intended use:
Marking microscope slides for purpose of identification.

Technical Specification :

Glass writing diamond tip
Handle made of wood, plastic, or Aluminium

Weight & volume:
Estimated weight: 0.011kg
Estimated volume: 1.00dm³

Packaging & labelling:
Unit presentation: Marker, diamond, set/3.



Laboratory Instrument India

