



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rack, test tube, 24 pos

Product Code :
LBNY-0007-2260002



Description :

Intended use:
To be used in the clinical laboratory setting.
To hold standard 18 mm Ø - 19 mm Ø tubes in an upright position.

Technical Specification :

Made of stainless steel
Provides positions to hold 24 tubes
Can hold standard 18 mm Ø - 19 mm Ø tubes
Stable base.

Weight & volume:
Estimated volume: 1.00dm³
Estimated weight: 0.09kg.



Laboratory Instrument India
