



Email : sales@laboratoryinstrumentindia.com

Product Name :

Rack, test tube, 24 pos

Product Code :

LBNY-0007-2260002



Description :

Intended use:

To be used in the clinical laboratory setting.

To hold standard 18 mm Ø - 19 mm Ø tubes in an upright position.

Technical Specification :

Made of stainless steel

Provides positions to hold 24 tubes

Can hold standard 18 mm Ø - 19 mm Ø tubes

Stable base.

Weight & volume:

Estimated volume: 1.00dm³

Estimated weight: 0.09kg.



Laboratory Instrument India

