



Email : sales@laboratoryinstrumentindia.com

Product Name :
Box, storage, 100 slides

Product Code :
LBNY-0007-2260003



Description :

Intended use:
Storage box for standard microscope slides with a size of 7.6 x 2.6 cm.

Technical Specification :

Flap cover has index card holder
Box made of resistant thermoplastic ABS
Space for 100 microscope slides
Box has hinged lid and a sturdy lock

Weight & volume:
Estimated volume: 0.14dm³
Estimated weight: 0.23kg.



Laboratory Instrument India