



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, sterile, screw cap, 2.0 ml, box/1000

Product Code :

LBNY-0007-2260004



Description :

Intended use:

Applicable for molecular/PCR & sterile microbiological works.

Various clinical laboratory usage conditions including; Sample cold storage, sample preparation/extractions with water bath, incubator, and centrifuge.

Technical Specification :

2.0ml capacity

Polypropylene material

Frosted side for Labeling

Skirted/Free standing base

Graduated

Knurled caps for firm grip when opening & closing

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Withstands centrifugation force up to 25,000 x g

Attached screw cap with O-rings to help prevent leaks.

Weight & volume:

Estimated volume: 12.6dm³

Estimated weight: 1.9kg.



Laboratory Instrument India