



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube, sterile, screw cap, 2.0 ml, box/1000

**Product Code :**

LBNY-0007-2260004

**Description :**

Intended use:

Applicable for molecular/PCR & sterile microbiological works.

Various clinical laboratory usage conditions including; Sample cold storage, sample preparation/extractions with water bath, incubator, and centrifuge.

**Technical Specification :**

2.0ml capacity

Polypropylene material

Frosted side for Labeling

Skirted/Free standing base

Graduated

Knurled caps for firm grip when opening & closing

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Withstands centrifugation force up to 25,000 x g

Attached screw cap with O-rings to help prevent leaks.

Weight & volume:

Estimated volume: 12.6dm<sup>3</sup>

Estimated weight: 1.9kg.



**Laboratory Instrument India**