



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube, sterile, screw cap, 5.0 ml, box/1000

**Product Code :**

LBNY-0007-2260005

**Description :**

Intended use:

Applicable for molecular/PCR & sterile microbiological works.

Various clinical laboratory usage conditions including; Sample cold storage, sample preparation/extractions with water bath, incubator, and centrifuge.

**Technical Specification :**

5.0ml capacity

Polypropylene material

Graduated

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Conical base

Screw cap

Withstands centrifugation force up to 25,000 x g

Weight & volume:

Estimated volume: 27.00dm<sup>3</sup>

Estimated weight: 5.00kg

---

## Laboratory Instrument India

