



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, sterile, screw cap, 5.0 ml, box/1000

Product Code :

LBNY-0007-2260005



Description :

Intended use:

Applicable for molecular/PCR & sterile microbiological works.

Various clinical laboratory usage conditions including; Sample cold storage, sample preparation/extractions with water bath, incubator, and centrifuge.

Technical Specification :

5.0ml capacity

Polypropylene material

Graduated

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Conical base

Screw cap

Withstands centrifugation force up to 25,000 x g

Weight & volume:

Estimated volume: 27.00dm³

Estimated weight: 5.00kg

Laboratory Instrument India

