



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, screw cap, 2.0 ml, box/1000

Product Code :

LBNY-0007-2260006



Description :

Intended use:

Various general clinical laboratory usage conditions including sample preparation, separation, centrifugation, and cold storage.

Technical Specification :

2.0ml capacity
Polypropylene material
Frosted side for Labeling
Graduated
Non-sterile
Skirted/Free standing base
Free of DNA, RNase, DNase & PCR inhibitors
Withstands centrifugation force up to 25,000 x g
Knurled caps for firm grip when opening & closing
Attached screw cap with O-rings to help prevent leaks

Weight & volume:

Estimated volume: 12.6dm³

Estimated weight: 1.9kg



Laboratory Instrument India