



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, screw cap, 5.0 ml, box/1000

Product Code :

LBNY-0007-2260007

**Description :**

Intended use:

Various general clinical laboratory usage conditions including; sample preparation, separation, centrifugation, and cold storage.

Technical Specification :

Polypropylene material

5.0ml capacity

Free of DNA, RNase, DNase & PCR inhibitors

Withstands centrifugation force up to 25,000 x g

Non-sterile

Conical base

Graduated

Screw cap

Weight & volume:

Estimated volume: 27.00dm³

Estimated weight: 5.00kg.

Laboratory Instrument India

