



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pipe Threading Machine Power

**Product Code :**  
LBNY-0004-820008



**Description :**

Pipe Threading Machine Power

**Technical Specification :**

Dies made of special pre-heat treated steel, ensuring long life and easy threading.  
Cutter & Reamer blade: most suitable shaped blade tip allows sharp cut and durability.  
Hammer chuck: Outer ring designed to grip easily and to grab pipes firmly with less power.  
The Jaw Body is wider giving greater load capacity and durability.  
Jaw insert improves the centripetal clamping feature.  
Spindle is supported by Sealed Ball Bearing. It ensures correct alignment and continued accuracy through lack of wear.  
Class F insulation grade on the motor allows continuous operation.  
The high capacity of Gear Pump operates even at low revolutions to assure precise threads.

**Specifications:**

Threading Capacity Water Pipe: 1/4"-3"  
Bolt: 3/8"2"  
Input Power: 1700W  
Output Power: 1140W  
Voltage: 240V,50Hz  
Spindle Speed: 33rpm

Accessory

---

Manual Die Head 1/2" to 2", 2Ea  
Manual Die Head 2-1/2" to 3", 2Ea  
BSPT Dies 1/2" to 3/4", 2Ea  
BSPT Dies 1"-1-1/2", 2Ea  
BSPT Dies 2" to 3", 2Ea



**Laboratory Instrument India**