




Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, screw cap, conic, 15/50ml,box/500

Product Code :

LBNY-0007-2260008



Description :

Intended use:

Various clinical laboratory usage conditions including Sample/reagent preparation, centrifugation, storage.

Technical Specification :

Polypropylene material

Screw cap, leak-proof

Sterile

Conical base, non-Skirted

Free of DNA, RNase, DNase & PCR inhibitors

Graduated

Supplied as a set of 2 different size tubes

1x Box of 250 tubes of 15ml

1x Box of 250 tubes of 50ml

Weight & volume:

Estimated volume: 150.00dm³

Estimated weight: 8.00kg.

Laboratory Instrument India

