



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette,transfer,sterile,3 ml,box/500

Product Code :

LBNY-0007-2260009

**Description :**

Intended use:

To transfer liquids in clinical laboratory settings applicable for molecular/PCR & sterile microbiological works.

Technical Specification :

Transfer pipette

Capacity 3.0 ml

Fine tip

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Embossed, minimal graduation 0.5 ml

With integrated flexible suction-bulb

Single-use disposable

Compatible with PCR applications & sterile microbiological applications

Individually peel-packed

Weight & volume:

Estimated volume: 16.5dm3

Estimated weight: 1.30kg.

Laboratory Instrument India

