



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipette, transfer, 3ml, box/500

Product Code :
LBNY-0007-22600010



Description :

Intended use:
To transfer liquids in general clinical laboratory applications.

Technical Specification :

Transfer pipette
Capacity 3.0 ml
Fine tip
Single-use disposable
Embossed, minimal graduation 0.5 ml
With integrated flexible suction-bulb

Weight & volume:
Estimated volume: 8.30dm³
Estimated weight: 0.75kg.



Laboratory Instrument India

