



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pipette, transfer, 3ml, box/500

**Product Code :**  
LBNY-0007-22600010



**Description :**

Intended use:  
To transfer liquids in general clinical laboratory applications.

**Technical Specification :**

Transfer pipette  
Capacity 3.0 ml  
Fine tip  
Single-use disposable  
Embossed, minimal graduation 0.5 ml  
With integrated flexible suction-bulb

Weight & volume:  
Estimated volume: 8.30dm3  
Estimated weight: 0.75kg.



**Laboratory Instrument India**

