




Email : sales@laboratoryinstrumentindia.com

Product Name :
Reservoir, reagent, 60ml, box/25

Product Code :
LBNY-0007-22600011



Description :

Intended use:
Reagent reservoir for short-time containment of reagent and chemicals during multichannel pipetting.

Technical Specification :

60ml capacity
Graduations on the inside wall
Polypropylene material, Re-usable, Autoclavable
V-bottomed trough
With pour-off spout
Wide base for stability
Compatible with the standard multi-channel pipettes (Up to 12 channels).

Weight & volume:
Estimated volume: 4.08dm³
Estimated weight: 0.75kg.



Laboratory Instrument India

