



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slide, microscope, frosted, box/100

Product Code :
LBNY-0007-22600012



Description :

Intended use:
Microscope slide for use with standard transmitted light microscope. For clinical laboratory use.

Technical Specification :

Material: Clear Microscopy optical grade glass
Clean, polished, free from nicks, scratches & fibres. Ready for use.
Sealed pack to prevent moisture contamination
Ground edges
Thickness: 1.0 to 1.2mm
Size: 76x26mm
Frosted end
Slides individually tissue paper interleaved.

Weight & volume:
Estimated volume: 0.36dm³
Estimated weight: 0.56kg.

Laboratory Instrument India

