



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slide, microscope, box/100

Product Code :
LBNY-0007-22600013



Description :

Intended use:
Microscope slide for use with standard transmitted light microscope. For clinical laboratory use.

Technical Specification :

Material: Clear Microscopy optical grade glass
Thickness: 1.0 to 1.2mm
Size: 76x26mm
Ground edges
Slides individually tissue paper interleaved
Sealed pack to prevent moisture contamination
Clean, polished, free from nicks, scratches & fibres. Ready for use.

Weight & volume:
Estimated volume: 0.36dm³
Estimated weight: 0.56kg.



Laboratory Instrument India

