



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cover glass, slides, box/100

Product Code :
LBNY-0007-22600014



Description :

Intended use:

Used as cover glass for microscope slides in standard transmitted light microscopy. For clinical laboratory use.

Technical Specification :

Material: Clear Microscopy optical grade glass

Clean, polished, free from nicks, scratches & fibres. Ready for use.

Thickness: 0.17mm

Size: 22x22mm

Sealed pack to prevent moisture contamination

Weight & volume:

Estimated volume: 4.00dm³

Estimated weight: 0.26kg.



Laboratory Instrument India

