



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tube,push cap,0.2 ml, ster,PCR,box/1000

**Product Code :**  
LBNY-0007-22600015



**Description :**

PCR tubes with push caps, 0.2ml, sterile, box/1000

Intended use:

PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

**Technical Specification :**

Single PCR tube, 0.2ml

Polypropylene material, thin wall with good thermal conductivity

Sterile, free of DNA/RNase/DNase and PCR inhibitors

With attached optical/flat cap

Compatible with standard format thermocycler blocks.

Weight & volume:

Estimated volume: 3.00dm<sup>3</sup>

Estimated weight: 0.23kg.



**Laboratory Instrument India**

