



Email : sales@laboratoryinstrumentindia.com

Product Name :

Rack,tubes,0.5/2.0/15.0ml,240 pos,set/3

Product Code :

LBNY-0007-22600016



Description :

Intended use:

To hold 0.5ml, 2.0ml microcentrifuge tubes/Cryovials, and 15ml tubes in an upright position. To be used in the clinical laboratory setting.

Technical Specification :

Made of rugged chemical resistant & autoclavable plastic

With alpha/numerical numbering of rack positions

Provides positions to hold 240 tubes

Holds 0.5ml, 2.0ml and 15ml tubes

Comes as a set of 3 racks.

Weight & volume:

Estimated volume: 2.80dm³

Estimated weight: 0.43kg.



Laboratory Instrument India

