



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Label,biohazard,adhesive2.5x2.5roll/1000

**Product Code :**  
LBNY-0007-22600017



**Description :**

Intended use:

Used in wide range of laboratory services from hospital laboratory to reference centre. Sticker for labelling biohazardous samples/waste materials.

**Technical Specification :**

Marked with English text "biohazard" and "black tri-sickle" logo (ref U+2623)

Strongly temper-proof self-adhesive

Colour: Red and black print and lining

Size: approximately 2.5 x 2.5 cm

Roll of 1000 labels.

Weight & volume:

Estimated volume: 4.00dm<sup>3</sup>

Estimated weight: 0.14kg.



**Laboratory Instrument India**

