



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Pressure Tester Power type

Product Code :
LBNY-0004-820009



Description :

Water Pressure Tester Power type

Technical Specification :

Water Pressure Tester (power type)

Specifications:

Electric Pressure Testing Pump

Tool Size: 250 * 530 * 410 mm

Standard Equipment: Pump + High Pressure Hose * Suction Pipe

Voltage: 230 VAC - 50 Hz

Pressure: 20 Bar

Flow Rate: 14-17 L/min

Tank: Epoxy Painted Steel

Connection Pipe: High Pressure 1 1/2" G Fem - 1500 Mm

Power: 750 W

Insulation Class; IP55

Gaskets: NBR

Laboratory Instrument India

