



Email : sales@laboratoryinstrumentindia.com

Product Name :

Box,stor0.5ml2.0ml5.0mltubes100pos,set14

Product Code :

LBNY-0007-22600018



Description :

Intended use:

Storage box for sample tubes, cryogenic vials and centrifuge tubes for clinical laboratory use.

Technical Specification :

Made of sturdy cardboard material with moisture repellent coating.

100-Cell divider for use with 1.5ml/2.0ml/5.0ml Cryovials/tubes.

Fits all standard stainless steel storage systems.

Withstands cold storage conditions including cryogenic freezers & liquid nitrogen.

With cover & cell-dividers.

Weight & volume:

Estimated volume: 32.00dm³

Estimated weight: 3.50kg.



Laboratory Instrument India

