



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Box,stor0.5ml2.0ml5.0mltubes100pos,set14

**Product Code :**

LBNY-0007-22600018



**Description :**

Intended use:

Storage box for sample tubes, cryogenic vials and centrifuge tubes for clinical laboratory use.

**Technical Specification :**

Made of sturdy cardboard material with moisture repellent coating.

100-Cell divider for use with 1.5ml/2.0ml/5.0ml Cryovials/tubes.

Fits all standard stainless steel storage systems.

Withstands cold storage conditions including cryogenic freezers & liquid nitrogen.

With cover & cell-dividers.

Weight & volume:

Estimated volume: 32.00dm<sup>3</sup>

Estimated weight: 3.50kg.



**Laboratory Instrument India**

