



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pipette, stand, 6 positions

**Product Code :**  
LBNY-0007-22600019



**Description :**

Intended use:  
Used to appropriately store multiple volumetric pipettes.

**Technical Specification :**

Made of chemical resistant material  
Bench-top stand, can hold 6 pipettes  
Sturdy construction for optimal stability  
Rotating as carousel  
Accommodates wide range of micropipettes

Weight & Volume:  
Estimated volume: 2.20dm<sup>3</sup>  
Estimated weight: 0.3kg.



**Laboratory Instrument India**

