



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette,tip,barri,ster,2-20ul,RACK-960

Product Code :

LBNY-0007-22600020



Description :

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

Technical Specification :

Variable Pipetting volumes 2-20ul
Sterile & free of DNA, DNase, RNase, and endotoxins
Smooth, Polypropylene material, prevents liquid retention
Racked primary packaging, 96 pipette tips per rack
Inert self-sealing barrier filter for prevention of aerosol contamination of pipette and sample
Universal fit, compatible with all standard micropipettes

Shelf life:
24-36 months

Weight & volume:
Estimated Volume: 0.01dm³
Estimated Weight: 6.01kg.

Laboratory Instrument India

