



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Pipette,tip,barri,ster,2-20ul,RACK-960

**Product Code :**

LBNY-0007-22600020

**Description :**

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

**Technical Specification :**

Variable Pipetting volumes 2-20ul

Sterile & free of DNA, DNase, RNase, and endotoxins

Smooth, Polypropylene material, prevents liquid retention

Racked primary packaging, 96 pipette tips per rack

Inert self-sealing barrier filter for prevention of aerosol contamination of pipette and sample

Universal fit, compatible with all standard micropipettes

Shelf life:

24-36 months

Weight & volume:

Estimated Volume: 0.01dm<sup>3</sup>

Estimated Weight: 6.01kg.

---

## Laboratory Instrument India

